

CAMERA & COMBO UNITS

Camera Combination System

REF: ES-12F1



This system is a combination of a high intensity Xenon light source, USB port for digit photo capture and download, a color CCD camera system, endoscope coupler, fiber-optic lighting cables and a 12 inch color LCD all enclosed in a compact portable design. The camera features an automatic shutter, single chip and is soakable. The Camera Combo Unit is universal voltage compatible for use anywhere in the world.

Micro Camera System

REF: EC-1S

The Micro Camera is a technically advanced single chip camera system. The automatic shutter ensures superior visualization. With a state of the art 1/2" chip that provides 500+ lines of resolution for the highest quality color and clarity. Camera comes complete with all cables and a 35mm low profile coupler.

