

ENDO-FOREHEAD

Endo-Forehead Punch

REF: EF-1010 Left Curved 12.5 cm

REF: EF-1011 Right Curved 12.5 cm

Allows atraumatic resection while minimizing the risk of damaging adjacent skin. The sharp, laterally-curved jaw is shielded by the lower biting surface when closed. Insulated shaft.



TUNGSTEN
CARBIDE

Endo-Forehead Grasping Forceps

REF: EF-2010 Left Curved 12.5 CM

REF: EF-2011 Right Curved 12.5 CM

Designed for manipulation and cauterization. Special jaws allow traction as well as grasping. Insulated shaft.



TUNGSTEN
CARBIDE

Nicanor Isse MD Endo-Forehead Elevator

REF: EF-3016-I

Designed by Nicanor Isse MD. 23 cm in length with finger rests for maximum control and comfort. Designed to make dissection easier for endo-facial procedures. Distal end has special angled, tapered tip.



STAINLESS
STEEL

ENDO-FOREHEAD

Endo-Forehead Elevator

REF: EF-3011 Half Curved

REF: EF-3012 Straight

REF: EF-3013 Quarter Curved

REF: EF-3014 Full Curved

This elevator is shaped to anatomical contours. The curvatures and tips are designed to give tactile feedback.



STAINLESS
STEEL

Endo-Forehead Nerve Elevator

REF: EF-3015

Tapered and sharp tipped to facilitate exposure.



STAINLESS
STEEL

Scalp Elevator

REF: EF-3010

Designed for soft tissue elevation. The flat tip with rounded edges permits atraumatic elevation while the curve permits easy access.



STAINLESS
STEEL

ENDO-FOREHEAD

Endo-Forehead Cannula

REF: EFC-0005 5mm with stop cock

REF: EFC-0004 4mm with stop cock

REF: EFC-0003 3mm with stop cock

REF: EFC-0002 2mm with stop cock

Specially designed for shielding the scope tip from surrounding tissues. Flat tip design for improved view.



STAINLESS
STEEL

“ISSE” Endo-Forehead Cannula

REF: EFC-0005-I 5mm “ISSE” Endo Cannula

REF: EFC-0004-I 4mm “ISSE” Endo Cannula

Designed by: Nicanor Isse MD.
A hooded distal end for improved view.



STAINLESS
STEEL