

FORCEPS

Adson 4 3/4 inches 4000 Jaw

REF: ESF-0500-TC

TUNGSTEN
CARBIDE

General purpose forceps
with wide gripping surface for
increased control.




Adson Dressing 4 3/4 inches Serrated

REF: ESF-0505

STAINLESS
STEEL

Stainless steel with serrated
jaws.




Adson Dressing 4 3/4 inches Cross Serrated

REF: ESF-0510

STAINLESS
STEEL

Stainless steel with cross
serrated jaws. Ideal for
removal of fine suture.



FORCEPS

Adson 1 x 2 Teeth 4 3/4 inches 4000 Jaw

REF: ESF-0515-TC

TUNGSTEN
CARBIDE

General purpose forcep with wide jaw surface and 1x2 teeth for increased control and gripping power.




Adson 1 x 2 Teeth 6.0 inches 4000 Jaw

REF: ESF-0520-TC

TUNGSTEN
CARBIDE

General purpose forcep with wide jaw surface and 1x2 teeth for increased control and gripping power.



Adson Tissue 1 x 2 Teeth 4 3/4 inches

REF: ESF-0525

STAINLESS
STEEL

Stainless steel with 1x2 teeth.



FORCEPS

Adson Tissue 1 x 2 Teeth 4 3/4 inches

REF: ESF-0530

STAINLESS
STEEL

Suture Tying Platform




Stainless steel with 1x2 teeth.

Brown Adson 4 3/4 inches 4000 Jaw

REF: ESF-0535-TC

TUNGSTEN
CARBIDE



Small multiple teeth eliminate traumatizing effect of the larger 1x2 teeth pattern.

Brown Adson 6.0 inches 4000 Jaw

REF: ESF-0540-TC

TUNGSTEN
CARBIDE



Small multiple teeth eliminate traumatizing effect of the larger 1x2 teeth pattern.

FORCEPS

Brown Adson 4 3/4 inches Brown Teeth

REF: ESF-0545

STAINLESS
STEEL



Stainless steel with multiple teeth. Ideal for atraumatic handling of tough tissue.

Gorney Adson 1x2 Teeth 4 3/4 inches

Diamond Dust

REF: ESF-0550-TC

TUNGSTEN
CARBIDE



Extremely fine 1x2 teeth with Diamond Dust tying platform for secure handling of delicate tissue and 6/0 to 8/0 suture.

Delicate 4 3/4 inches Diamond Dust

REF: ESF-0555-TC

TUNGSTEN
CARBIDE



Designed for accurate tissue and needle handling. Fine teeth and Diamond Dust tying platform provides secure grip of fine needles and tissue. Ergonomic design and non-slip handle for comfort and control.

FORCEPS

Castroviejo 4.0 inches Diamond Dust

REF: ESF-0560-TC



TUNGSTEN
CARBIDE

.3mm tip

Extremely delicate tissue forceps with fine teeth and Diamond Dust tying platform to enhance handling of delicate tissue and suture.



Castroviejo 4.0 inches Diamond Dust

REF: ESF-0565-TC



TUNGSTEN
CARBIDE

.6mm tip

Extremely delicate tissue forceps with fine teeth and Diamond Dust tying platform to enhance handling of delicate tissue and suture.



Bayonet 7 3/4 inches 7000 Jaw

REF: ESF-0615-TC

TUNGSTEN
CARBIDE

Rounded handle bayonet forceps designed for maximum visualization and fingertip control. Features a carbide insert for secure needle handling and tissue grasping.



FORCEPS

Julian Suture 4 3/8 inches Smooth Jaw

REF: ESF-0575

STAINLESS
STEEL

Excellent for removal of delicate sutures due to smooth jaw on a raised platform.



Bishop Harmon 3 1/2 inches 1x2 Teeth

REF: ESF-0580

STAINLESS
STEEL

Light action forceps ideal for grasping of delicate tissue.



Gillies Dissecting 6.0 inches Serrated

REF: ESF-0585

STAINLESS
STEEL

Fine jaws and delicate 1x2 teeth provide secure grip.



FORCEPS

Cushing 7.0 inches 4000 Jaw

REF: ESF-0590-TC

TUNGSTEN
CARBIDE



Ideal for soft tissue work and general suturing.

Potts Smith 8.0 inches 4000 Jaw

REF: ESF-0595-TC

TUNGSTEN
CARBIDE



Useful for grasping delicate tissue and suture.

Potts Smith 9.0 inches 4000 Jaw

REF: ESF-0600-TC

TUNGSTEN
CARBIDE



Useful for grasping delicate tissue and suture.

FORCEPS

Gerald 7.0 inches 7000 Jaw

REF: ESF-0605-TC

TUNGSTEN
CARBIDE

Delicate tapered 1mm tip for working with fine tissues and suture.



DeBakey 8.0 inches 1.5mm Jaw

REF: ESF-0610

STAINLESS
STEEL

DeBakey serration for the atraumatic handling of delicate tissue.

