

RETRACTORS

Breast Retractor Fiber-Optic & Suction

REF: BR-5538 Serrated blade 90 mm x 24 mm

REF: BR-5539 Serrated blade 90 mm x 30 mm

REF: BR-5540 Serrated blade 115 mm x 24 mm

REF: BR-5541 Serrated blade 150 mm x 36 mm

REF: BR-5542 Serrated blade 190 mm x 30 mm

REF: BR-5543 Smooth blade 90 mm x 24 mm

REF: BR-5544 Smooth blade 90 mm x 30 mm

REF: BR-5545 Smooth blade 115 mm x 24 mm

REF: BR-5546 Smooth blade 150 mm x 36 mm

REF: BR-5547 Smooth blade 190 mm x 30 mm

STAINLESS
STEEL



Retractors designed to enhance visualization and control during a wide range of procedures. Suction port to provide for smoke evacuation if needed. Serrated blades prevent tissue slippage.

**ALSO AVAILABLE WITH OUT FIBER-OPTICS,
CALL FOR MORE DETAILS AND ORDER NUMBERS**

RETRACTORS

Double Ended Breast Retractor

REF: BR-5610

STAINLESS
STEEL

Versatile double-ended retractor for multiple applications. Ergonomically designed handle maximizes comfort while using retractor for long periods of time.



Economic Breast Retractor Serrated

REF: ER-4380

STAINLESS
STEEL

Dual sided economical breast retractor with teeth on each end to prevent tissue slippage. Made of 400 series surgical stainless steel.




“S” Retractor

REF: ER-4370

STAINLESS
STEEL

General “S” shaped retractor. Made from Hardened 400 series surgical stainless steel.



RETRACTORS

Biggs Face Lift

REF: ER-4600

STAINLESS
STEEL

Superior retraction at the distal end, while allowing visualization. The curve removes the surgeons hand from the operating field.



Freeman Flap

REF: ER-4590 Width 1"

STAINLESS
STEEL

Tissue is held securely but gently in the bend of the retractor. Tip is slightly curved downward for elevation and visualization.




Benicke

REF: ER-5640 19.0 cm Sharp

REF: ER-5641 19.0 cm Blunt

STAINLESS
STEEL



A break through in double retractors. A uniquely designed "push", "pull" hook at one end allows for use in multiple procedures.



RETRACTORS

Aufricht

REF: ER-4620

REF: ER-4610 With Fiber Optics

Designed with a custom curved and atraumatic stabilizing teeth at the tip for stability.



STAINLESS
STEEL

Mini Langen Beck

REF: ER-4730 10 x 6 mm

Thin blade general retractor. Tip slightly angled to prevent slippage.



STAINLESS
STEEL

Langen Beck

REF: ER-4740 30 x 11 mm

General retractor. Tip slightly angled to prevent slippage.



STAINLESS
STEEL

RETRACTORS

DesMarres Lid

REF: ER-4920 5 1/2" x 10mm

REF: ER-4930 5 1/2" x 10mm Insulated

REF: ER-4940 5 1/2" x 12mm

REF: ER-4950 5 1/2" x 12mm Insulated

REF: ER-4960 5 1/2" x 14mm

REF: ER-4970 5 1/2" x 14mm Insulated

REF: ER-4980 5 1/2" x 16mm

REF: ER-4990 5 1/2" x 16mm Insulated

STAINLESS
STEEL



Popular general purpose retractors used on delicate tissues.

Deaver Retractor

REF: ER-4800 1 inch x 9 inches

STAINLESS
STEEL



Popular general purpose retractors used on delicate tissues.

Deaver Retractor

REF: ER-4810 1 inch x 12 inches

REF: ER-4820 1 inch x 13 inches

REF: ER-4830 1 1/2 inch x 12 inches

REF: ER-4840 2 inch x 12 inches

REF: ER-4850 3 inch x 12 inches

STAINLESS
STEEL



Popular general purpose retractors used on delicate tissues.