

# MARKERS

## McKissock Keyhole Pattern

REF: ER-4710 50mm Pattern


REF: ER-4720 45mm Pattern

REF: ER-4630 42mm Pattern

REF: ER-4640 38mm Pattern

Accurately marks amount of tissue for excision. Limbs are semi-malleable for fine adjustments.

STAINLESS  
STEEL



## Freeman Areola Marker

REF: ER-4650 36mm Marker

REF: ER-4660 34mm Marker

REF: ER-4670 50mm Marker

REF: ER-4680 45mm Marker

REF: ER-4690 42mm Marker

REF: ER-4700 38mm Marker

Facilitates marking of the areola prior to incision.

STAINLESS  
STEEL



## Four Ring Marker

REF: EM-4540

Allows for choice of four different diameters; 40mm, 44mm, 48mm and 52mm. Must be used with marking ink.



# MARKERS

## D'Assumpcao F/L

REF: EM-2300

STAINLESS  
STEEL



Facilitates the precise marking of dissected flaps for excision. Eliminates the guess work.

## D'Assumpcao F/L with Ratchet

REF: EM-2310

STAINLESS  
STEEL



Facilitates the precise marking of dissected flaps for excision. Eliminates the guess work.

## Umbilicas Marker

REF: ER-5642 290 mm 180mm

STAINLESS  
STEEL



Allow for the precise, error free marking of the umbilicas.