

# MEDICAL MONITORS

## 17 Inch Flat Screen Medical Grade

REF: EX-17WS



17 inch LCD with a 4:3 aspect ratio. This unit is water sealed with an anti-glare glass front. Capable of 16.7 million colors and a resolution of 1280 x 1024. The 17 inch monitor has a viewing angle of 170 degrees and a response time of 15/10 ms. This unit is cart or wall mountable and has video inputs for VGA, Composite or S-Video. The monitor is compatible with NTSC, PAL or SECAM.

---

## 19 Inch Flat Screen Medical Grade

REF: EX-19WS



19 inch LCD with a 4:3 aspect ratio. This unit is water sealed with an anti-glare glass front. Capable of 16.7 million colors and a resolution of 1280 x 1024. The 19 inch monitor has a viewing angle of 170 degrees and a response time of 15/10 ms. This unit is cart or wall mountable and has video inputs for VGA, Composite or S-Video. The monitor is compatible with NTSC, PAL or SECAM.