

Modular Laparoscopic

Handle Options

Selecting a modular laparoscopic instrument.

Begin by selecting the handle type that you desire from one of the pictures. Below the picture you will notice a description and part number, the part number of the handle is the first portion of your instruments order number. Next, select the jaw style you desire from one of the pictures. Below the picture you will notice a description and part number again, the part number of the jaw is the next portion of your instruments order number.

(e.g. : K225.000-K221.5340.1 = Rotating metal handle w/ no ratchet, 5mm Metzenbaum scissors straight).

Standard shaft length is 33cm, plus length of jaw. Custom lengths available upon request.



Metal Handle Isolated
Angled With HF Rotate

Part Number: K220.0000



Plastic Handle With HF
Rotate

Part Number: K221.0000 (Ratchet)
Part Number: K221.0000.O (No Ratchet)



Metal Handle Isolated
Angled Rotate
(No Ratchet)

Part Number: K225.0000



Metal Handle Isolated
Angled Rotate
(Ratchet)

Part Number: K226.0000

Modular Laparoscopic

Shaft Tubes Isolated 33 cm



5 MM Diameter Shaft Tube
Isolated


Part Number: K221.5300

10 MM Diameter Shaft Tube
Isolated

Part Number: K221.1300


Modular Laparoscopic

Assembly



Modular Laparoscopic

Disassembly



Modular Laparoscopic

Scissors: 5 mm Diameter




Matzenbaum Straight

Part Number: K221.5340.I




Matzenbaum Curved

Part Number: K221.5341.I



Micro Dissection Straight

Part Number: K222.5301.I




Micro Dissection Curved

Part Number: K222.5302.I


Modular Laparoscopic

Scissors: 5 mm Diameter



Hook


Part Number: K221.5312.I



Peritoneal Blunt


Part Number: K221.5314.I

Punches: 5 mm Diameter



Biopsy Without Thorn

Part Number: K221.5310.I




Biopsy With Thorn

Part Number: K221.5311.I

Modular Laparoscopic

Forceps: 5 mm Diameter



Tumor Grasping Forceps Pointed

Part Number: K221.5318.I



Biopsy Without Thorn

Part Number: K221.5301.I



Biopsy With Thorn

Part Number: K221.5302.I




Atraumatic

Part Number: K221.5305.I


Modular Laparoscopic

Forceps: 5 mm Diameter




Atraumatic With Reservoir

Part Number: K221.5309.I




Dolphin

Part Number: K221.5315.I



Conical With Reservoir

Part Number: K221.5316.I




Dissection Fenestrated Slim

Part Number: K221.5307.I

Modular Laparoscopic

Forceps: 5 mm Diameter



Concave Smooth Blunt

Part Number: K221.5306.I



Sharp Cross Teeth Jaw 18MM

Part Number: K221.5343.I



Sharp Cross Teeth Jaw 14MM

Part Number: K221.5344.I




Dissection 4:2 Teeth Double Action

Part Number: K221.5338.I


Modular Laparoscopic

Forceps: 5 mm Diameter




Dissection 4:2 Teeth Single Action

Part Number: K221.5352.I




Dissection 2:1 Teeth

Part Number: K221.5327.I



Dissection Distal With Row of Teeth

Part Number: K221.5339.I




Maryland Dissector Curved

Part Number: K221.5321.I


Modular Laparoscopic

Forceps: 5 mm Diameter



Extraction Retrograde Serrated

Part Number: K221.5335.I




Universal

Part Number: K221.5317.I



Claw Grasper 2:3 Teeth

Part Number: K221.5308.I




Alligator Jaw

Part Number: K221.5322.I


Modular Laparoscopic

Forceps: 5 mm Diameter



Dissection Fenestrated

Part Number: K221.5329.I



Dorsey

Part Number: K221.5349.I



Babcock

Part Number: K221.5319.I




Allis With Debakey Teeth

Part Number: K221.5342.I


Modular Laparoscopic

Forceps: 5 mm Diameter



Davis & Geak Wavy Tooth W/O Distal Teeth

Part Number: K221.5323.I



Davis & Geak Wavy Tooth W/ Distal Teeth

Part Number: K221.5324.I



Johann Duck Bill

Part Number: K221.5350.I



Mixer

Part Number: K221.5332.I

Modular Laparoscopic

Needle Holders: (Carbide) 5




Tungsten Carbide Needle Holder Straight

Part Number: K221.5345.I



Tungsten Carbide Needle Holder Left Curve

Part Number: K221.5346.I



Tungsten Carbide Needle Holder Right Curve

Part Number: K221.5347.I


Modular Laparoscopic

Forceps: 10 mm Diameter



Claw Forceps 2:3 Teeth

Part Number: K221.1008.I



Extraction Retrograde Serrated

Part Number: K221.1035.I



Maryland Dissector Curved

Part Number: K221.1001.I



Mixer

Part Number: K221.1032.I

Modular Laparoscopic

Forceps: 10 mm Diameter

**Image Not
Available**

Dorsey

Part Number: K221.1048.I




Babcock Grasping Fenestrated

Part Number: K221.1019.I



Allis With Debaquey Teeth

Part Number: K221.1042.I



Babcock With Cross Teeth

Part Number: K221.1043.I

**Image Not
Available**

Tumor Grasping Forceps

Part Number: K221.1037.I